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The enclosed BEGO-GOLD card is part of these processing instructions!

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Alloy characteristics

Biocertificate	<input checked="" type="checkbox"/>	Melting interval [°C]	1140–1055
Type	hard (3)	Casting temperature [°C]	1320
BEGO-GOLD colour code	yellow 4	Preheating temperature [°C]	850
Density [g/cm³]	18.0	CTE 20–600 °C (25–500 °C)	
Average grain size [µm]	20	[10°K¹]	14.6 (14.4)
Vickers hardness (HV 5)	* 130/175/190	Heat treatment	500 °C 15 min.
Elongation limit (Rp 0.2) [MPa]	* 225/430/480	Soft annealing	750 °C 10 min.
Ductile yield (A5) [%]	* 13/9/8	(then quenching in water at 20 °C)	
Modulus of elasticity approx. [MPa]	92,000	* soft / after firing / hardened	

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Processing instructions

Modelling:

- Minimum metal thickness (after grinding): for ceramic veneering 0.4 mm, for acrylic veneering with retention pearls 0.3 mm.
- Connecting parts between the pontics should be as thick and high as possible (at least 3.5 mm high and 2.5 mm wide). Provide palatal interdental garlands as reinforcement for wide span bridges.
- Spruing of single crowns: provide casting reservoir.

Investing: Use phosphate-bonded crown and bridge investment materials (e.g. Auroplus®, Bellavest®).

Casting and finishing: General: Do not overheat alloy. Use only clean crucibles, one crucible per alloy. Recommendation: to enable an exact identification of each case cast new metal only

- In case of re-casting: only re-cast identical alloys. Blast old material. Add at least 50 % of new material.
 - Use ceramic or graphite crucible. Ceramic crucible: sprinkle some Auromelt HF melting powder on the casting ingots.
 - Continue to heat after complete melting of the casting pieces:
- | | |
|-----------------------|-------------------|
| Flame melting: | 3–5 seconds |
| HF induction heating: | approx. 3 seconds |
| Resistance heating: | 2–3 minutes. |

Use fine carbide or BEGO sintered diamond milling tools for finishing.

Ceramic: Use veneering ceramics in accordance with DIN EN ISO 9693 with firing temperatures of up to approx. 960 °C (e.g. Carat, Biodent, Duceram, IPS-Classic, Omega, Omega 900, VMK 95). Always follow the ceramic manufacturer's instructions!

- Support frames carefully (use refractory wool, additional wires or make individual firing supports of BegoForm®).
- Set the heating rate to a maximum of 55 °C per minute!

Oxide firing:

- Always blast the surface to be veneered (Korox® 110, 2 – max. 3 bar) and clean the frame thoroughly (steam clean or boil in aqua dest.).
- Maintain oxide firing at 950 °C without vacuum for 10 minutes (Omega 900: 900 °C).
- Acid treatment of oxide with Aurocid (60 °C, 10–15 minutes). Clean thoroughly (steam clean or boil in aqua dest.).

Ceramic firing: Long-term cooling is recommended.

Soldering:

- Support object in a soldering block of Bellatherm®. Prepare a gap of max. 0.2 mm with parallel walls.
- Soldering before firing with the flame (1030 °C): PontoStar® G-Solder (Order no. 61045) and Minoxid or Fluxsol.
- Soldering after firing in furnace (810 °C): BEGO-Gold-Solder I (Order no. 61017) and Minoxid. Allow to cool normally.
- Acid-treat residual flux in Aurocid (60 °C, approx. 1 minute). Clean thoroughly (steam clean or boil in aqua dest.).

Laser welding: Filling material: PontoStar® G wire, dia. 0.4 mm (Order no. 61150) or cast one yourself (Acrylic profile dia. 0.45 mm, Order no. 53005).

Warranty: Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products undergo constant further development and are therefore subject to modification regarding design and composition.

Secondary effects: Like e.g. allergies to contents of the alloy or electrochemically based reactions may very rarely occur.

Reciprocal actions: In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Reactions: In case of known incompatibilities and allergies to contents of the alloy.